

SCBI for STEAM with 8 hour PH reaction



Quickly monitoring with SCBI

The monitoring of sterilization condition at preferred places with the help of Self Contained Biological Indicators (or SCBI) is no longer practice of every day in European Hospitals and mainly used during validation of the sterilizers, whereas outside Europe the SCBI is widely used for daily routine checks of sterilizers in Hospitals.

REF: 107.010.0100
REF: 120.001.0001

SCBI for STEAM, pop. 1×10^6
Dual Incubator for SCBI

Unit: 100 pieces
Unit: 001 piece

Characteristics:

- ◆ Biological Indicators conform to ISO 11138
- ◆ SCBI STEAM: *Geobacillus Stearothermophilus* ATCC 7953 1×10^6 spores per vial
- ◆ Color change: Purple to Yellow (STEAM)
- ◆ Digital Incubator (220V) for 60°C / 8 hrs (STEAM)