

## Steam Integrator - Type 5 - Self-adhesive



Integrators can be used in addition to the batch control Helix which will allow you to release load. Especially when OR staff wants to be re-affirmed about the sterilization conditions inside trays or containers.

REF: 105.203.0500 Steam Integrator 134°C-3,5 min Unit: 500 pieces 121°C-15 min

Steam Integrators are giving a clear result which easily can be interpreted by the end-user by comparing it to the reference color printed on the integrator. No misunderstanding possible.

## Characteristics:

- ◆ Steam Integrator conform to ISO 11140 Type 5
- ◆ Clear area to write details like dates, department, operator
- ◆ Shelf life: 5 years after production
- ◆ Packed in an easy dispensing carton
- ◆ Clear color transition from blue into pink
- ◆ Lacquer covered indicator
- ♦ Non-Toxic / Lead Free





BIC : TGBATRIS Branch No.: 437
IBAN: TR36 0006 2000 4370 0009 0896 01
IBAN: TR09 0006 2000 4370 0009 0896 02



Tel : +90-212-613 80 54

Fax : +90-212-613 80 55

Email: info@spmedikal.com

Web: www.spmedikal.com

Bank: Garanti Bank - Istanbul - TURKEY

(EUR) Account No.: 9089601

(US\$) Account No.: 9089602

VAT: Bayrampaşa - 781 046 3595

BIC: TGBATRIS Branch No.: 437

IBAN: TR36 0006 2000 4370 0009 0896 01

IBAN: TR09 0006 2000 4370 0009 0896 02