

Steam Indicator - Type 4



The multi-parameter indicator is checking time and temperature of the sterilization cycle at the exact place where the indicator is placed. A clear Violet to Green color change is indicating the correct steam sterilization cycle parameters have been met.

REF: 104.003.0500 Multi-parameter indicator

134°C-3,5 min / 121°C-15 min Unit: 500 pieces

The Sterintech™ Indicators are intended to be used during steam sterilization to check critical parameters at the place where they are placed.

Characteristics:

- ♦ Shelf life: 5 years after production
- ◆ Packed in an easy dispensing carton
- ◆ Steam indicator conform to ISO 11140 :2014 Type 4
- ♦ Clear color transition into green
- ◆ Non-Toxic / Lead Free

Other types and lay-outs available upon request.

