

SCBI for EO - 12 hour PH reaction





Quickly monitoring with SCBI

The development towards quicker read-outs of SCBI is still ongoing , most of the solutions require a Quick Read-out Incubator which is expensive to acquire. The new Sterintech SCBI EO is giving a PH reaction with a normal dry bed incubator in 12 hours time only. An improvement compared to our previous model which was needing 48 hours of incubation.

REF: 107.030.0100 SCBI for EO, pop. 1x10E6 Unit: 100 pieces
REF: 120.001.0001 Dual Incubator for SCBI Unit: 001 piece

Characteristics:

♦ Biological Indicator conform to ISO 11138 and ISO 11138-2

♦ SCBI EO: Bacillus Atrophaeus ATCC 7953

→ Population: 1x10E6 spores per vial
 → Color change: Orange to bright Yellow





BIC : TGBATRIS Branch No.: 437
IBAN: TR36 0006 2000 4370 0009 0896 01
IBAN: TR09 0006 2000 4370 0009 0896 02