

## SCBI for STEAM, Plasma & Ethylene-Oxide



## Quickly monitoring ..... with SCBI

The monitoring of sterilization condition at preferred places with the help of Self Contained Biological Indicators (or SCBI) is no longer practice of every day in European Hospitals and mainly used during validation of the sterilizers, whereas outside Europe the SCBI is widely used for daily routine checks of sterilizers in Hospitals.

 REF: 107.001.0100
 SCBI for STEAM, pop. 1x10E6
 Unit: 100 pieces

 REF: 107.002.0100
 SCBI for Plasma, pop. 1x10E6
 Unit: 100 pieces

 REF: 107.003.0100
 SCBI for EO, pop. 1x10E6
 Unit: 100 pieces

 REF: 120.001.0001
 Dual Incubator for SCBI
 Unit: 001 pieces

## Characteristics:

- → Biological Indicators conform to ISO 11138
- ◆ SCBI EO: Bacillus Atrophaeus ATCC 9372 1x10E6 spores per vial
- ◆ SCBI STEAM/Plasma: Geobacillus Stearothermophilus ATCC 7953 1x10E6 spores per vial

Bank: Garanti Bank - Istanbul - TURKEY

(EUR) Account No.: 9089601

(US\$) Account No.: 9089602

VAT: Bayrampaşa - 781 046 3595

◆ Color change: Green to yellow (EO) and purple to yellow (STEAM)

Tel: +90-212-613 80 54

Fax : +90-212-613 80 55

Email: info@spmedikal.com

Web: www.spmedikal.com

◆ Digital Incubator (220V) for 37°C / 48 hrs (EO) and 60°C / 24 hrs (STEAM & Plasma)



IBAN: TR36 0006 2000 4370 0009 0896 01

IBAN: TR09 0006 2000 4370 0009 0896 02

BIC: TGBATRIS Branch No.: 437